
ACCU-CHEK® FastClix

Lancing device
Autopiqueur
Dispositivo de punción
Lancetador
採血筆
قلم الوخز

Instructions for use

The Accu-Chek FastClix lancing device is a device for simple, virtually pain-free collection of capillary blood from the fingertips. The lancing device has 11 different penetration depths so that you can adjust the penetration depth to the texture of your skin. With the release button, you can prime the lancing device and trigger a lancet in a single step. When a lancet is triggered, it penetrates the skin. You insert a lancet drum with 6 lancets into the lancing device. You can obtain blood using a new, sterile lancet 6 times before the lancet drum needs changing.

These instructions for use feature the following 3 symbols:

- This symbol indicates a **possible risk of injury or of damage to your health**.
- This symbol draws attention to actions that could result in **damage to the lancing device**.
- This symbol draws your attention to **important information**.

• The Accu-Chek FastClix lancing device is **intended for personal use only**! It is only allowed to be used to obtain blood from one person.

There is a risk of infections being transmitted if the lancing device is used by other people, even by family members, or if healthcare professionals use this lancing device to obtain blood from different people. This lancing device is therefore not suitable for professional use in healthcare facilities and institutions.

- Keep the lancing device with all its components away from children under the age of 3 years. There is a risk of suffocation if small parts (e.g. caps or similar objects) are swallowed.

Only insert Accu-Chek FastClix lancet drums into the Accu-Chek FastClix lancing device. Other lancets or lancet drums may impair the operation of the lancing device or damage it permanently.

Preparation


- Wash your hands with warm water and soap before you obtain blood (see 1). Dry your hands thoroughly. This reduces contamination of the puncture site and stimulates blood flow.

Inserting the lancet drum

You have to insert a lancet drum to be able to obtain blood with the lancing device.

- Remove the cap from the lancing device (see 2).
- Insert a new lancet drum – white end first – into the lancing device (see 3).

The lancet drum is properly inserted when it is completely inside the lancing device, and only the projecting edge protrudes from the lancing device at the front (see 3).

You must not insert the lancet drum into the lancing device and simultaneously press the release button or hold the lancing device such as a table top. This could release a lancet and inadvertently cause injury.

- Place the cap back on again. The cut-out section at the edge of the cap must be around the correspondingly shaped curvature with the penetration depth indicator (silver square) (see 4).

The cap is on correctly when it is pushed back as far as it will go, and you hear and feel it click into place.

The cap can only be correctly placed on up to the stop position when the lancet drum is completely inside the lancing device. Do not use any force to place the cap on as this might damage the cap and the lancing device. Correct the position of the lancet drum instead.

Each lancet drum contains 6 lancets. The first lancet is now ready to be used. The window showing the number of available lancets shows 6; there are now 6 lancets available (see 5).

- Do not remove the lancet drum before all 6 lancets have been used. The lancet drum cannot be reused once it has been removed from the lancing device.

Setting the penetration depth

You can set 11 penetration depths (0.5 to 5.5). The silver square (penetration depth indicator) on the lancing device indicates the current penetration depth setting. The larger the number, the greater the penetration depth.


If you have no experience with this lancing device, we recommend setting a low penetration depth, such as penetration depth 2.

- Rotate the cap until the desired penetration depth lines up with the indicator (see 6).

Obtaining a blood drop

In general, blood can be obtained from any finger. Certain fingers can be unsuitable if, for example, a skin or fingernail infection is present. We recommend obtaining capillary blood from the sides of the fingertips as these areas are the least sensitive to pain.

- Always use a new lancet when you obtain blood. In doing this, you reduce the risk of infection, and obtaining blood remains virtually pain free.
- Only use the lancing device when the cap is attached. Without the cap, the lancet penetrates too deeply and the puncture can become painful.

• In individual cases, the lancet may protrude from the cap after obtaining blood. There is a risk of injury or infection if you or another person touches the protruding lancet. If the lancet protrudes from the cap, **carefully** set the lancing device aside. Keep the lancing device out of reach of other people and contact customer support.

- Press the lancing device firmly against the selected puncture site (see 7).
- Press the release button all the way down (see 8).

The lancet is released and it penetrates the skin.

- Massage the finger in the direction of the fingertip to encourage a blood drop to form. The amount of blood that emerges depends on the penetration depth and the pressure used to hold the lancing device against the skin.

If not enough blood emerges, apply more pressure on the lancing device when you next obtain blood. If that is not sufficient, increase the penetration depth progressively in addition.

- If too much blood emerges, decrease the penetration depth.
- Test your blood glucose immediately after you obtained blood, according to the instructions for use for your blood glucose meter.

We recommend protecting the puncture site from contamination afterwards.

Loading a new lancet

If you have already used the lancing device to obtain blood, you will need a new lancet the next time you obtain blood.

- Push the lever at the side of the lancing device as far as it will go in direction A and then back in direction B (see 9).

The next lancet is now ready to be used.

The window next to the lever shows the new number of available lancets, which has decreased by one, for example, 5 instead of 6.


• Proceed as in the previous step each time you need to change the lancet.

Each time a new lancet is loaded, the number of available lancets shown in the window decreases by one, for example, from 4 to 3. When you have loaded the sixth and last lancet, the window shows 1.

Replacing a used lancet drum

When you have used the sixth and last lancet, you need to replace the used lancet drum with a new one.

- Remove the cap from the lancing device (see 2).
- Remove the lancet drum from the lancing device (see 10).

When a lancet drum has been used, a red stripe is visible on the white cylinder.

- Insert a new lancet drum (see section "Inserting the lancet drum").

A used lancet drum has a built-in lock that stops you from using it again.

Used lancet drums can be disposed of in household waste* if no other regulations apply in your country. *In Australia subject to local council regulations.

Storing the lancing device

Do not store the lancing device at extreme temperatures, for example, in a hot car. This can impair the operation of the lancing device.

Cleaning and disinfecting the lancing device and cap

To prevent infections from being transmitted, you must clean and disinfect the lancing device and cap regularly.

- At least once a week.
 - When there is blood on the lancing device or cap.
 - Always before someone else handles the lancing device, for example, to assist you.
- You need lint-free cloths in which you can wrap the whole lancing device, water and 70 % isopropanol. Before you disinfect the lancing device and the cap, you first have to clean them to remove blood and other dirt.

- To avoid damaging the lancing device and the cap:
- Do not dip the lancing device and the cap in the fluids.
- Use only water for cleaning and only 70 % isopropanol for disinfecting.


- Remove the cap from the lancing device (see 2).
- Remove the lancet drum from the lancing device (see 10).

Cleaning

- Moisten a cloth with water.
- Wipe the lancing device and the cap thoroughly; the lancing device from the outside only (see 11) and the cap from the outside and the inside (see 12).
- Dry the lancing device and the cap with a dry cloth.

Disinfecting

- Dampen a cloth with 70 % isopropanol.
- Wipe the lancing device and the cap thoroughly for 2 minutes each; the lancing device from the outside only and the cap from the outside and the inside.

- Place the cap back on again.
- Wrap the whole lancing device in the cloth.
- Leave the lancing device wrapped for 8 minutes.
- Unwrap the lancing device and let it dry in the air.

Alternative sites for obtaining blood

You can obtain capillary blood for blood glucose testing not only from the fingertips, but also from other sites (palm, forearm, upper arm, calf or thigh), so-called alternative site testing (abbreviated as AST). To be able to obtain capillary blood from other sites, you have to replace the regular cap by a special cap, the Accu-Chek FastClix AST cap. Before you decide to use blood from alternative sites, consult your healthcare professional.

Ordering information

If you want to order an Accu-Chek FastClix lancing device, Accu-Chek FastClix lancet drums or an Accu-Chek FastClix AST cap, contact customer support (see section Customer Support).

Last update

2019-04

Customer Support


Roche Diabetes Care Australia Pty. Limited
2 Julius Avenue, North Ryde
NSW 2113, Australia
ABN 69 602 140 278
Accu-Chek Enquiry Line: 1800 251 816
www.accu-check.com.au

Roche Diabetes Care
division of Hoffmann-La Roche Limited
201 Armand-Frappier Blvd.
Laval (Québec)
H7 V 4A2
Canada
Accu-Chek Customer Care: 1-800-363-7949
www.accu-check.ca

Hong Kong
Enquiry hotline: +852-2485 7512 (office hours)
www.accu-check.com.hk

Malaysia
Customer Support and Service Centre:
Toll-Free Line: 1-800-88-1313
www.accu-check.com.my

Roche Diabetes Care South Africa (Pty) Ltd.
Hertford Office Park, 90 Bekker Road
Vorna Valley, 1686
South Africa
Tel: +27 (11) 504 4600
www.accu-check.co.za

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
Thailand
Customer Service Line: +66 (0) 2791 2222
www.accu-check.co.th

	Caution, refer to safety-related notes in the instructions for use accompanying this product.
	Use only once
	Use by
	Sterilized using irradiation
	Manufacturer
	Catalogue number
	Batch code
	Global Trade Item Number
	This product fulfils the requirements of the European Directive 93/42/EEC on medical devices.

ACCUCHEK and FASTCLIX are trademarks of Roche.

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© 2020 Roche Diabetes Care